



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

RAD51AP1-RAD51 associated protein 1 Gene

Reference: PTXBC016330

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PIR51; RAD51-associated protein 1; RAD51-interacting protein; RAD51 associated protein 1		

For research use only.

GENE DESCRIPTION	RAD51 associated protein 1
GENE ID	10635
ORIGIN SPECIES	Human
ACCESSIONS	BC016330
VECTOR	pDONR223
NCBI SYMBOL	RAD51AP1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgttgcggcctgtgagacataagaacacagctcaattactcacagttgtaccactctgacagtgatgatgtttgtttctgcaactctgaccttaaacagaacaaatccagaacagcaccaaaggagttaaaaacaagataaaccaaacctactctgaaacaatctccggaagaagaaatccag
tacaagaacaaaccccttaaaaaagagactccttggaagacttttagttactgcacagctgcagctgcttgcacaaaagatgctcttagatgacagctctacacagagacttagaagttgctgactgctttatcagtgaggagaaactccacactcagctcaccatatttgcgcgaactctcagataa
agacttgaaaacactgcgcatgtgtataaataagaacaattgaataactctctccataactcttaattgcgtgtgacagctgattatgattgttgataagattactctggaaagatgattgttggtgtgtcagaggaagaaacacagactctaaagctcagcagcagcagaagaattcttctg
gaaggcagtgatggtgtgatgtcttaatgacactgacaccagacttgcacctgtgtgaagattctgaggatgattctgtttgtgagatgtgaggataatgacgaagactctctatgagaagaagaataagtagtaagaataataaaagaagaagtgagggttaaaatccccagtagaaaaa
gagaagaatactaaactcaaatgtaatgtgttggtgactctgggtgagctctgcacgtgcgctgcacaaatcagaactcagctccttgcctcaaaaagagttctctgctctcagatcagactgaaaacacttagaataatcagcagctctcagctgacgtcagctgacgtcagcaagaacacttaaatgggtgccaccagc
ggcatctgtgtgcagagaaatgcacagccacactggtgtatgtctgtgaaatgcctcaacaaatgctccgccttgcgtctgcagatgacacagagctgaaactctgtcatccaaatgccactcagctga

More info: Please contact us for any question: contact@proteogenix.fr