

## CCNL2-cyclin L2 Gene

Reference: PTXBC016333

|                |                                                                                                                                        |                     |                |
|----------------|----------------------------------------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                                                    | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                                                   | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | ANIA-6B; CCNM; CCNS; HCLA-ISO; HLA-ISO; PCEE; SB138; cyclin-L2; cyclin M; cyclin S; paneth cell-enhanced expression protein; cyclin L2 |                     |                |

**For research use only.**

---

|                         |                                                                                                   |
|-------------------------|---------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | cyclin L2                                                                                         |
| <b>GENE ID</b>          | 81669                                                                                             |
| <b>ORIGIN SPECIES</b>   | Human                                                                                             |
| <b>ACCESSIONS</b>       | BC016333                                                                                          |
| <b>VECTOR</b>           | pDONR223                                                                                          |
| <b>NCBI SYMBOL</b>      | CCNL2                                                                                             |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                    |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon )

atggcggcgggcgggcgggcggtgctgctgagggcggcagctcccgcgcgagcgggcgggcggtgctggtggcgccctcagggtcgagggggtgctgctgctggcgagggcgtgactccgggtgctcatcacttgagaaactgcctcctgctgacgacaa  
gctccgtttcacgccgtccatgtcgagcggcctcgacagcgacacagagaccgacctccgctggtgggtgctgagctcctcaggcgggcggtatcctgctccgctgcccaggtggccatggctaccggcgaggtgtgtccagcgggtctttataccaagtccctcgtggaagcac  
tccatggagcatgtgtcaatggcctgtgtccacgtggcttccaaatagagaagggcccaagacgcatacgggacgtcatcaatgtgtttaccgcttcgacagctgagagacaaaaagaaagccgtgctctactactggatcaagattatgtaattaaagaacaaattataaaggcg  
gaaagacgagttctcaaaaggtgggtttctgctccatgtgaagcatcctcataagataatcgttatgtacctcaggtgttagaggtgagcgtaaccaacacctggtccagacctcatggtagccttgagggtgaagtga