

ZNF509-zinc finger protein 509 Gene

Reference: PTXBC016477

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ZNF509; zinc finger and BTB domain-containing protein 49; zinc finger protein 509; zinc finger and BTB domain containing 49		

For research use only.

GENE DESCRIPTION	zinc finger protein 509
GENE ID	166793
ORIGIN SPECIES	Human
ACCESSIONS	BC016477
VECTOR	pDONR223
NCBI SYMBOL	ZNF509
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcaaggaggttagtaaaacacagaattcggcacacggggagcgggcttacagctgctctgctcgggaaatgttttggggatcaggtgacctcgcaggcatgtccgcactcacactggggagaagccgtacacatgtgagatctgtaacaagtgtttacccgctctcggtgctccggcgccacaagaagatgcactgcaaaagctggtgacgagagccagatgtgctggaggagctcagccaagccatcgagacctccgacctcgagaatctcagagctcagactctttctcccaagacacgtctgtgacgctgatgccagtgtcgggttaaactccctgtccaccagtgaaaattctgtgtagaattgatagccactctggcggtctctattgttaagtacggtccatgatccaactcatggagtttagtgaccaggagaagctgagtttgatcctgggtaaacttgccaagccccagatgcagcagacacagcctcaggcctatgcttactcggatgtggacccccagccggtggcgaaccactgcaggccgatggcatggccatgatcgttctctgctgtgcttggacaaccacggcggtgacccccctgggcagtcgagcatcttcaccacttataggaaactcagagggtcagttttctccagcatgactctctggggctagcgtgaagacgctgcagaatgaaaacgagtttagaccagtga