

RAB37-RAB37, member RAS oncogene family Gene

Reference: PTXBC016615

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	RAB37, member RAS oncogene family; RAB37, member of RAS oncogene family; ras-related protein Rab-37		

For research use only.

GENE DESCRIPTION	RAB37, member RAS oncogene family
GENE ID	326624
ORIGIN SPECIES	Human
ACCESSIONS	BC016615
VECTOR	pDONR223
NCBI SYMBOL	RAB37
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgacgggcacgcaggcgcgttgccacccgggatggcagggccccgagcgcctccccgcctgcagtcgagctacgacctacgggcaaggatgcttctgggagacacaggcgtcggcaaaacatgtttcctgatccaattcaagacgggacctctgtccggaaccttc
atagccaccgtcggcatagactcaggacaagggtggtgactgtggtgagcgtgagagtgaaagctgcagatctgggacaccgtgggcagggaacggttcgaagcgtcacccatgcttattacagagatgctcaggcctgtctctgtgtatgacatcaccacaattcttctgacaa
catcaggccctggctcactgagattcatgagtatgccagaggacgtggtgatcatgctgctaggcaacaaggcgatagcagcgaagagtgatccgttcgaagacggagagaccttggccagggtacggtgttccctctggagaccagcgccaagactggcatgaat
gtggagttagccttctggccatcgccaaggaaactgaataccgggcccggcatcaggcggatgagccagcttcagatccgagactatgtagtcccagaagaagcgctccagctgctgctcctcatgtga