

SSX5-synovial sarcoma, X breakpoint 5 Gene

Reference: PTXBC016640

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	protein SSX5; synovial sarcoma, X breakpoint 5; SSX family member 5		
For research use only.			

GENE DESCRIPTION	synovial sarcoma, X breakpoint 5
GENE ID	6758
ORIGIN SPECIES	Human
ACCESSIONS	BC016640
VECTOR	pDONR223
NCBI SYMBOL	SSX5
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaacggagacgatgcctttgtacggagacctaggggtgtctcaataaccacagaagatgcaaaagcatccctggagacaagtctgtgaccgtggaatacatttggatgaatcagtcctgttggaaaggtgggaagagagccagccagcagcattaaagctctactgtgtggcaggg
gagaagctagggccttcgatgataatgccaaatacttcttgagaagagtgaggaaagatgaaagcctcggagaaaatcatctatgtgtatatgaagagaaagtatgagccatgactaaactaggtttcaaggccacctccacctttcatgcgtaataaacgggtcgagacttcagg
ggaatgatttgataatgacctaacctgggaatcaggttgaacatcctcagatgacttcggcaggtccagggaatcttccgaagatcacgcccagagaagccagcagaggaaggaaatgattcaaaaggagtgccagaagcatctggcccacagaacaatgggaacagctgcgc
ccctcaggaactaaatacctctgagaaggttaacaagacatctggacccaaaaggggaaacatgcctggacccacagagtgctgtgagagaaagcaactggtgatttatgaagatcagcgacctcaggaagatgacgagtaa