

FOSL1-FOS-like antigen 1 Gene

Reference: PTXBC016648

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	FRA; FRA1; fos-related antigen 1; FOS like 1, AP-1 transcription factor subunit; FOS-like antigen-1; FOS like 1, AP-1 transcription factor subunit		

For research use only.

GENE DESCRIPTION	FOS-like antigen 1
GENE ID	8061
ORIGIN SPECIES	Human
ACCESSIONS	BC016648
VECTOR	pDONR223
NCBI SYMBOL	FOSL1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgttcgagacttggggaacccggcccgagctcgggaacggcggtacggcgcccgccgagcccgccgagcagcagcagcagcagaagttccactggtgccaagcatcaacacatgagtgagtcaggagctgagtgatggtacagcctattcc
tggggcccgagcagttaccccgagcctctgacctaccctcagtagacccccacaccccgccagggagtcacccggccctggggccgctccaggggtacgtcgaaggcctgtgaacagatcagcccgagggaagaggagcgcgcggaacaag
ctggctgcccgaagtgcagggaaccggaaggaactgaccgacttctgcaggcggagactgacaaactggaagatgagaactctgggctgcaggagagattgaggagctgcagaagcagaaggagcgccctagagctggtgctggaagcccaccgacccatctgcaaaatc
ccggaaggagccaaggggggacacagcagtagcagtgccaccagcagccaccagcccgcccgccctgtgaccttgccttccctccagggcctgtgcttgaacctgagcactgcacacccacactatgaccacacccctcctaactcttccccagcctggt
cttcacctacccagcactcctgagccttggctcagctcgcgaagtagcagcagcagcgagaccatccttgaccccttggctctccaacccctcctcgttctga