



Reference: PTXBC016652

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	BMX non-receptor tyrosine kinase; Etk/Bmx cytosolic tyrosine kinase; cytoplasmic tyrosine-protein kinase BMX; ETK; PSCTK2; PSCTK3; BTK-like on X chromosome; NTK38 tyrosine kinase; bone marrow tyrosine kinase gene in chromosome X protein; epithelial and en		

For research use only.

GENE DESCRIPTION	BMX non-receptor tyrosine kinase
GENE ID	660
ORIGIN SPECIES	Human
ACCESSIONS	BC016652
VECTOR	pDONR223
NCBI SYMBOL	BMX
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr