

TSLP-thymic stromal lymphopoietin Gene

Reference: PTXBC016720

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	thymic stromal lymphopoietin		

For research use only.

GENE DESCRIPTION	thymic stromal lymphopoietin
GENE ID	85480
ORIGIN SPECIES	Human
ACCESSIONS	BC016720
VECTOR	pDONR223
NCBI SYMBOL	TSLP
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaaaactaaggctgccttagctatctgtgtcccaggtatcggaaactcagataaatgctactcaggcgaatgaagaagaggagaaaaggaaagtcacaccaaataatgtctggaacaagtgtcacaattacaaggattgtggcgtcgttcaatcgaccttactgaacaacagtaa