

CAPZA3-capping protein (actin filament) muscle Z-line, alpha 3 Gene

Reference: PTXBC016745

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CAPPA3; Gsg3; HEL-S-86; F-actin-capping protein subunit alpha-3; CP-alpha-3; CapZ alpha-3; F-actin capping protein alpha-3 subunit; capping protein (actin filament) muscle Z-line, alpha 3; epididymis secretory protein Li 86; germ cell-specific protein 3;		

For research use only.

GENE DESCRIPTION	capping protein (actin filament) muscle Z-line, alpha 3
GENE ID	93661
ORIGIN SPECIES	Human
ACCESSIONS	BC016745
VECTOR	pDONR223
NCBI SYMBOL	CAPZA3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgacacttagcgtgctgagcaggaaggacaaggaaagagtaattcgacactgtttattacagggccctccagggaattgttaaatgctttgatgatctctgtctgcttaccgtgatgaaaaacttatgcaccaccaaggtgagtggtgcaggccaccaactgcaaaaaattctgtacc
actctgcatcgatggaatccagtactctgtctaccacaatgtatggcgactaccgatttttgaccatcaaaagcaaaactttcttcaaatattacgtcttcaaaatcagctgaaagacatccaaagtcagtgatcattcagaatgaggcagaatacctgagagttgttctgtgctgctta
aaactgtatgtgaatgaccactatccaaaaggaaatgcaacatgctgagaaaaactgtcaaaagttaaggagcttgaatgcttgaatgaagatcacaactatgaacaggagagtgctggaacgactttggaatctaaatggaatttccaagttaaccatttcaaccagtaacggg
aagaatattgtgcaagctcacttctcagggtgtgcaacctcatattgaaatatccaaaggacctgaaagaaagcttggaatagtttaaccaagctcaactggctctaagtttgcaggctgtggaagagcaagagaacaatttcaagctgcagctctggaagaattacagggtatccaat
gaagccctgagaaaaattctacgaaggatcttccagtgcaccgcactcttattgactggcagagatactctgacttgaatctggtgatgtatcctaaataggatgtcatttattcaagaagtggtgtgcaactggataataa