

C14orf147-chromosome 14 open reading frame 147 Gene

Reference: PTXBC016805

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C14orf147; SSSPTA; serine palmitoyltransferase small subunit A; small subunit of serine palmitoyltransferase A		

For research use only.

GENE DESCRIPTION	chromosome 14 open reading frame 147
GENE ID	171546
ORIGIN SPECIES	Human
ACCESSIONS	BC016805
VECTOR	pDONR223
NCBI SYMBOL	C14orf147
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgctggcgcggcctggaagcagatgtcctgttctactaccagtagctgtgtcacggcgtctacatgctggagccctgggagcggagcgtgttcaattccatgctgtttccattgtggggatggcactatacacaggatagcttctcatgccccagcacatcatggcgatattg
cactactttgaaatcgtacaatga