



Reference: PTXBC016853

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	3-hydroxy fatty acids omega-hydroxylase CYP4F11; phylloquinone omega-hydroxylase CYP4F11; CYP1VF11; cytochrome P450 4F11; cytochrome P450, family 4, subfamily F, polypeptide 11; cytochrome P450, subfamily IVF, polypeptide 11; cytochrome P450 family 4 subf		

GENE DESCRIPTION	cytochrome P450, family 4, subfamily F, polypeptide 11
GENE ID	57834
ORIGIN SPECIES	Human
ACCESSIONS	BC016853
VECTOR	pDONR223
NCBI SYMBOL	CYP4F11
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

Page 1/1