

TRIM31-tripartite motif-containing 31 Gene

Reference: PTXBC016866

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	E3 ubiquitin-protein ligase TRIM31; C6orf13; HCG1; HCGI; RNF; tripartite motif-containing protein 31; tripartite motif containing 31		

For research use only.

GENE DESCRIPTION	tripartite motif-containing 31
GENE ID	11074
ORIGIN SPECIES	Human
ACCESSIONS	BC016866
VECTOR	pDONR223
NCBI SYMBOL	TRIM31
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaataaaatgacatgaagagctggggctgtttacagaaaaatacataaaatgaacaaacctcagagcccggtcatcttctgcaggcggcagaactacatcggggccacaaatcaccactctcagcccatccactccctgtttcgggctcgtctgctgggaaagtcacttttccagtatgtctctggcctcttatgatgagatttctggtcaaggagcagctctcaggatcacgaagacatttgacgttgcgctgtccgaggagctccatgcggcactgagtgagtgctgacagcgatccgggcttggtttgtgaggttccttcaagctaa