



Reference: PTXBC016906

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	SEC12; prolactin regulatory element-binding protein; mammalian guanine nucleotide exchange factor mSec12; prolactin regulatory element binding		

For research use only.

GENE DESCRIPTION	prolactin regulatory element binding
GENE ID	10113
ORIGIN SPECIES	Human
ACCESSIONS	BC016906
VECTOR	pDONR223
NCBI SYMBOL	PREB
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]