

PRR5-proline rich 5 (renal) Gene

Reference: PTXBC016921

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	FLJ20185k; PP610; PROTOR-1; PROTOR1; proline-rich protein 5; Rho GTPase activating protein 8; proline rich 5 (renal); protein observed with Rictor-1; proline rich 5		

For research use only.

GENE DESCRIPTION	proline rich 5 (renal)
GENE ID	55615
ORIGIN SPECIES	Human
ACCESSIONS	BC016921
VECTOR	pDONR223
NCBI SYMBOL	PRR5
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggccccatgccacaggacagaagctgctggactcactggcagagacctgggacttcttcaagtacgtgctgccccatgctcagggccatcttaccgggtgcaggccaaggagccatcggtgcgccagctggccctgctgcacttccggaatgccatcacccctcagtgtaagctagaggatgcgtggccggcccatgccgcgtgccccctgccatcgtgcagatgctgctgctgcaggggtacatgagtcaggggcgtgactgaggactacgtgcctggagacgtgggtccagaaggtggtgctgccatactgggcacctacggcctccactccagcgagggcccttacccttctgcatcctggaagagcgcctctccgccgtccctcgcggggacgtgctggccaagaacctgtggtgctgctccaagagctacaacacgcctctgctgaaccccggtgcaggagcacgagggcgagggcgccggcgccggcggtaccagcatccgcaggcactctgtgctggagatgacgtcctgccccagcctcagggtctctccgacccggccagggcccccacgggacgttcaggtctccccggcgccccactcaggccctgccccagcagactgtacccccagccccagccccctgagcaggggttgatccacccgcagctccctgccccgtccagccccgagaacctggtggaccagatcctggagtcctggactcggattctgaaggattttcattgacttggccggggcggtcctgcatgtccgacttgagggccttgggggcgagagtgctgtgta