



Antibodies and Proteins

ProteoGenix SAS
Espace Européen de l'Entreprise
15 rue de La Haye
67300 Schiltigheim, France
(0)3 90 20 54 70 - Fax +33 (0)3 88
contact@proteogenix.fr
www.ProteoGenix.science

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88
contact@proteogenix.fr
www.ProteoGenix.science

TLL1-tolloid-like 1 Gene

Reference: PTXBC016922

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ASD6; TLL; tolloid-like protein 1; tolloid like 1		
For research use only.			

GENE DESCRIPTION	tolloid-like 1
GENE ID	7092
ORIGIN SPECIES	Human
ACCESSIONS	BC016922
VECTOR	pDONR223
NCBI SYMBOL	TLL1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggggtgggaacgctttcccgaggatgctgtgtggctgtggctctggggtattgttcacggggagctatgggtctgcgctggcctgattatgattacactttgatgggaacgaagaggataaaacagagactatgattacaaggaccctgttaaagcgcgtgtatttggggcg
atattgccttagatgtaggaacacttaataacttcaaaatagatggaaactgacttaccagacacccctttggaaacctggacataccacagggttggaagactgctatgccaagaagcagggggccctctcaacactctatagacagatagaagaattgttcgttctggag
caaaaacacacagttgaaagaaatgactctcaaaattctcaggccaaatgagaanaaactgagttccagacgctctacatcaagaagggaataatggctctggaggctgtttctctatgtatagaagaagaaactctgcgacccaagaagcagatgtcaagcagccatgag
cactctggaaaagacacatgctgtgacttttcaataaagaagatgafgaaagactatattgacttccattatgagatctgttgatgtctgctctatctgtatgtctgcggagaaagatcgagacgaactctatctgcgaagactgtgataatttgggtgttctgag
tagctttctgcgatgaacacacagacatcgatgagataaccgatgactatcaagaagaataactcagccaggctcaagatgaactatttctgaagtgtgagctctgagagctgaagaaataaactcattggaagaagaatgattgttcagactatgctcactatccggaagcaactctcaagg
ggatgtttctgataccattctcccctcctgtagtaaatggctagactcttccaaatgttgcacgacacccctctaaagcaaaaggagatgctgcacaggcaagaagcgtgtatagtgctccagcatgttggaagaaactctcaagaatccaatggcaactttctctcaggatttcccaatggc
tacccttctccagctccatcgaaagattctctctaccccgacgaagatgatttttctacctttga