

AKAP7-A kinase (PRKA) anchor protein 7 Gene

Reference: PTXBC016927

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	AKAP15; AKAP18; A-kinase anchor protein 7 isoform gamma; A kinase (PRKA) anchor protein 7; A-kinase anchor protein 18 kDa; A-kinase anchor protein 7 isoforms alpha and beta; A-kinase anchor protein 9 kDa; AKAP 18; AKAP-7 isoform gamma; AKAP-7 isoforms alp		

For research use only.

GENE DESCRIPTION	A kinase (PRKA) anchor protein 7
GENE ID	9465
ORIGIN SPECIES	Human
ACCESSIONS	BC016927
VECTOR	pDONR223
NCBI SYMBOL	AKAP7
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggccagcttctgcttcttctcaagagatgaaggaaaaatcagtgaaaagaacggagggagcccgatgacgctgaactagtaaggctcagtaagagcgtgttgagaaacgcggtgctcaaggctgtccagcagtagtctggaggaaacacagaataaaaaaagccggggg
aggggagctctgtgaaaaccgaagcagctgatcagaatggcaatgacaatgagaacaacaggaaatga