

FLJ35767-FLJ35767 protein Gene

Reference: PTXBC016939

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	testis-expressed sequence 19 protein; testis expressed 19		

For research use only.

GENE DESCRIPTION	FLJ35767 protein
GENE ID	400629
ORIGIN SPECIES	Human
ACCESSIONS	BC016939
VECTOR	pDONR223
NCBI SYMBOL	FLJ35767
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtgccctccggtcagcatcggtatgaggaagagggcatgtctacgtctacgtctcctggatgtatcagcttcaacatggagatcagctaacatttgcctcacctgcttcaaggctgccttctagactttaagacttgcctgagtcagaggactgggaagaagacaactgggacctg
agctgatggagcacactgagcagagtcagagcaggggggtcctcaggatggagctgagctggggcgagagccagcagcctgtgcaggggggctctgagcagtgaggccctggcagcagccccagaagggttgaagatgcaggtctgacccccacttt
gtccccactgaactatggcctcaggagctgtgccccctgggctgggcttgaggatgctgactggaccagggtcttcctggagatttgaggagcttctacctgctcacactggccaagcttcttctctcatag