

NLRP1-NLR family, pyrin domain containing 1 Gene

Reference: PTXBC016965

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CARD7; CIDED; CLR17.1; DEFCAP; DEFCAP-L/S; NAC; NALP1; PP1044; SLEV1; VAMAS1; NACHT, LRR and PYD domains-containing protein 1; NACHT, LRR and PYD containing protein 1; NACHT, leucine rich repeat and PYD (pyrin domain) containing 1; NACHT, leucine rich rep		

For research use only.

GENE DESCRIPTION	NLR family, pyrin domain containing 1
GENE ID	22861
ORIGIN SPECIES	Human
ACCESSIONS	BC016965
VECTOR	pDONR223
NCBI SYMBOL	NLRP1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgccctggacccctaccactctgtgacgtgggggcacgcgcaggatcccatcttctgtgttggggagctcagagtgcgccaatcttgaatcttaaggatgagccagaccagaccgcggcctctagagaggtccggcaggagggtcggcgcctggccgggggtg
ggccggaggccctgtgatgctgcatcgcccccaggaggagccagctgtgcccagagttggcggccgagagagagacaagagcgcgcagcaggcgaagctggaggcgggacgcgacttctgtcgtgcccggaggagtcgagactgtacccggaggagctgtctcaccag
gagaccacgtcctggaagtgtccgggacgcgcgggacgtgtggctgcagaccccgccggcagcagggccagagctggcgactcctga