

NECAB2-N-terminal EF-hand calcium binding protein 2 Gene

Reference: PTXBC016979

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	EFCBP2; stip-2; N-terminal EF-hand calcium-binding protein 2; EF-hand calcium-binding protein 2; neuronal calcium binding 2; neuronal calcium-binding protein 2; synaptotagmin-interacting protein 2; N-terminal EF-hand calcium binding protein 2		

For research use only.

GENE DESCRIPTION	N-terminal EF-hand calcium binding protein 2
GENE ID	54550
ORIGIN SPECIES	Human
ACCESSIONS	BC016979
VECTOR	pDONR223
NCBI SYMBOL	NECAB2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggtataccaagaaggtatatgaggggtgggagcaacgtggaccagttgtgacccgttctcctgaaggagacggccaatcagatccagtcgtgctgagctcagtgagagtgccggtggagggccatcgaggaaacagaccagccagctccgacagaaccacatcaaacccagc
cacagcgcggcacagacctggtgtggaagccccactccgcctctgcccccaaccacaagctcatggtatggaacaaggcaagaccctccatctgccacggaggtgcaaaaggaaaggggtctggaagccagatcagccgcttgccagagctgattgggaggtggagagcaa
agcactgtggttcgacctgcagcagcgccgtgcatgaagatggcaccaacatgcacctgcagctggtccggcaggagatggcgtgtgccccgagcaactgagcgagtttctggactctctgcgcagatctgcggggaccactggcgtgaggaactgctccacatcactgccg
tgaggctctcagatggcttaccctttgtcatctatgagttctgggagacagaggagcgtggaagaggcactgcagagccccctgtgtaaggcgttccggcacgtcaaggtggacacactgagccagcctgaggccctctccagatcttggtccagctgcttggcaggtgggac
gggactga