



Reference: PTXBC017023

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	MBD1; MDP; RDP; dipeptidase 1; dehydropeptidase-I; hRDP; microsomal dipeptidase; renal dipeptidase; testicular tissue protein Li 57; dipeptidase 1 (renal)		

For research use only.

GENE DESCRIPTION	dipeptidase 1 (renal)
GENE ID	1800
ORIGIN SPECIES	Human
ACCESSIONS	BC017023
VECTOR	pDONR223
NCBI SYMBOL	DPEP1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]