



Reference: PTXBC017033

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	squalene monooxygenase		
For research use only.			

GENE DESCRIPTION	squalene epoxidase
GENE ID	6713
ORIGIN SPECIES	Human
ACCESSIONS	BC017033
VECTOR	pDONR223
NCBI SYMBOL	SQLE
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

Page 1/1