



Reference: PTXBC017045

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	GPBB; glycogen phosphorylase, brain form; phosphorylase, glycogen; brain; glycogen phosphorylase B		

For research use only.

GENE DESCRIPTION	phosphorylase, glycogen; brain
GENE ID	5834
ORIGIN SPECIES	Human
ACCESSIONS	BC017045
VECTOR	pDONR223
NCBI SYMBOL	PYGB
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]