



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

GPR175-G protein-coupled receptor 175 Gene

Reference: PTXBC017058

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	integral membrane protein GPR175; GPR175; TMEM227; TPRA40; transmembrane protein adipocyte-associated 1; G protein-coupled receptor 175; seven transmembrane domain orphan receptor; transmembrane domain protein regulated in adipocytes 40 kDa; transmembrane		

For research use only.

GENE DESCRIPTION	G protein-coupled receptor 175
GENE ID	131601
ORIGIN SPECIES	Human
ACCESSIONS	BC017058
VECTOR	pDONR223
NCBI SYMBOL	GPR175
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr