



Reference: PTXBC017061

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	transmembrane GTPase MFN2; CMT2A; CMT2A2; CMT2A2A; CMT2A2B; CPRP1; HMSN6A; HSG; MARF; mitofusin-2; hyperplasia suppressor; mitochondrial assembly regulatory factor; mitofusin 2		

For research use only.

GENE DESCRIPTION	mitofusin 2
GENE ID	9927
ORIGIN SPECIES	Human
ACCESSIONS	BC017061
VECTOR	pDONR223
NCBI SYMBOL	MFN2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]