

C3orf52-chromosome 3 open reading frame 52 Gene

Reference: PTXBC017064

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	TTMP; TPA-induced transmembrane protein; TPA induced trans-membrane protein; novel endoplasmic reticulum transmembrane protein; chromosome 3 open reading frame 52		

For research use only.

GENE DESCRIPTION	chromosome 3 open reading frame 52
GENE ID	79669
ORIGIN SPECIES	Human
ACCESSIONS	BC017064
VECTOR	pDONR223
NCBI SYMBOL	C3orf52
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgatcatcaccctccatttctaggtgtcattacagtgatcatcataggttatgtcttgcagtaacttatgttgatgaagatgaaatgaaatactgaattatcatcaacaaaacattctcatcatgctgaagattccagaggagtggttgcgaagagggaattgcctcacctgctaccga
aaggctcacagatgtgacagtacatcgccctctctgagtcgttatttacttcagttgaaatagtggaacttcagtggtgaaaafgccacagtaacgtatgacctgcaattgggttccatcagatgatgaaaatttatgaagtatatgatgagtgaggagttggtgctgggcattttgctacagga
ttccgtgatcagaataacctggttggagagctctggggcttgatccaacatccctcttgcctatgaatga