

ORMDL3-ORM1-like 3 (*S. cerevisiae*) Gene

Reference: PTXBC017087

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ORM1-like protein 3; ORMDL sphingolipid biosynthesis regulator 3		
For research use only.			

GENE DESCRIPTION	ORM1-like 3 (<i>S. cerevisiae</i>)
GENE ID	94103
ORIGIN SPECIES	Human
ACCESSIONS	BC017087
VECTOR	pDONR223
NCBI SYMBOL	ORMDL3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaatgtgggcacagcgacagcgaggtgaacccaacacgagggtgatgaacagccgtggcatctgctctctacgtgctggccatcgggtctctccacatcgtgctgagcagatcccggtttgtgagtgctcctgtcgtctggaccctcaccaacctcattcacaacatggcatgtat
atcttctgcacacggtgaaggggacacctttgagacccggaccagggcaaggcgaggtgctaaccactgggagcagatggattatgggtccagttcagggcctctcggaagttcttgaccatcacaccatcgtgctgtacttctcaccagcttctacactaagtacgaccagat
ccattttgtgtcaacacggtgtccctgatgagcgtgttatccccaagctgcccagctccacggagtcgggatttttggatcaataagtactga