

POMZP3-POM (POM121 homolog, rat) and ZP3 fusion Gene

Reference: PTXBC017101

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	POM-ZP3; POM121; POM121 and ZP3 fusion protein; POM (POM121 homolog, rat) and ZP3 fusion; POM (POM121 rat homolog) and ZP3 fusion; POM-ZP3 fusion protein; POM121/ZP3 fusion protein; POM121 and ZP3 fusion		

For research use only.

GENE DESCRIPTION	POM (POM121 homolog, rat) and ZP3 fusion
GENE ID	22932
ORIGIN SPECIES	Human
ACCESSIONS	BC017101
VECTOR	pDONR223
NCBI SYMBOL	POMZP3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggtgtgtagccagtgactctgaggatcgccctcctgacagaagatttcggttctgcgataccagagcagataatcagctcaacactgtctcaccatcaagtaatgccagaccatgtgcaaggagactgtactgagtgccctcaagagaagaagaagaaaggacagtgg
aggagaagaccaatattccttgatggccaggaaaataaagaagctgtctgtcagcgtctcactgatgcctcttctgcatcacaagtctcctgacccggccagatacactccagttcacagtgatctctccacttgctaatgactccagaacatgatatacatcacctgccactga
agggtcaccttagctgagcaggaccagatgaactcaacaaggcctgttcttcagcaagccttccaacagctggttccagtggaaggccggctgacatctgtcaatgctgtacaagggtgactgtggcactccaagccattccaggaggcagcctcgtgctgtgagccagtgtgtcca
cgtctgttcccgtaacccgagcatgtgacagaagaagcagatgtcaccgtggggccactgactctcttgacagcagtggtga