

POLR2I-polymerase (RNA) II (DNA directed) polypeptide I, 14.5kDa Gene

Reference: PTXBC017112

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	RPB9; hRPB14.5; DNA-directed RNA polymerase II subunit RPB9; DNA-directed RNA polymerase II 14.5 kDa polypeptide; DNA-directed RNA polymerase II subunit I; RNA polymerase II 14.5 kDa subunit; RNA polymerase II subunit B9; RPB14.5; polymerase (RNA) II (DNA		

For research use only.

GENE DESCRIPTION	polymerase (RNA) II (DNA directed) polypeptide I, 14.5kDa
GENE ID	5438
ORIGIN SPECIES	Human
ACCESSIONS	BC017112
VECTOR	pDONR223
NCBI SYMBOL	POLR2I
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggagcccgacgggacttacgagccgggttcgtgggtattcgttctgcccaggaatgaacaacatgctgtaccccaaggagagacaaggaacgcattctgctctacgcgtgccggaactgtgattaccagcaggagccgacacagctgcatctatgtcaacaagatcacgca
cgaagtggagcgaactgaccagattatcgccgacgtgtccaggaacccacgttgcgcgggacggaggaccacccgtgccaaaagtgcggccacaaggaggctgtgtttctccagtacacagtgccgggcccaggagccatgcgcctttactacgtgtgcacagccccacactg
cggccaccgctggaccgagtga