

## GPC4-glypican 4 Gene

Reference: PTXBC017166

|                               |                                                                                     |                     |                |
|-------------------------------|-------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:                 | DNA                                                                                 | Storage Condition:  | -20°C          |
| Size/Quantity:                | 2 ug                                                                                | Delivery Condition: | +4°C, Blue Ice |
| Synonyms                      | K-glypican; glypican-4; dJ900E8.1 (glypican 4); glypican proteoglycan 4; glypican 4 |                     |                |
| <b>For research use only.</b> |                                                                                     |                     |                |

|                         |                                                                                                  |
|-------------------------|--------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | glypican 4                                                                                       |
| <b>GENE ID</b>          | 2239                                                                                             |
| <b>ORIGIN SPECIES</b>   | Human                                                                                            |
| <b>ACCESSIONS</b>       | BC017166                                                                                         |
| <b>VECTOR</b>           | pDONR223                                                                                         |
| <b>NCBI SYMBOL</b>      | GPC4                                                                                             |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                   |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon)

atggcacgggttcggtgcccgcgtctctgcacctggcagtgctcagcgcgcgctgctggcgtccgagctcaagtgcaaaagtgtcgcgaagtgcgacgtctttacgtgtccaaaggcttcaacaagaacgatgccccctccacgagatcaacggatgatcatttgaagatctgtcc  
ccagggttctacgtgctctcaagagatggaggagaagtacagcctgcaaaagtaagatgattcaaaagtgtggtcagcgaacagtgcaatcatttgaagctgtctttgcttcacgttacaagaagtttgatgaattcttcaaaagaactacttgaaaatgcagagaatccctgaatgatat  
gtttgtgaagacatattggccatttatacatgcaaaattctgagctatttaagatctctctgtagagttgaaacgttactacgtgtgtgggaaatgtgaacctgggaagaaatgctaaatgacttctgggtcgcctcctggagcggatgttccgctgtggaactccagtgaccatttacagatgagt  
atctggaatgtgtgagcaagtatacggagcagctgaagcccttcggagatgtccctcgcaaatggaagctccaggttactcgtgtttgttagcagcccgactttcgtcaaggcttagcgtgtgctgggagatgtctgtgagcaaggtctcctgtggaaccccccacagccagtgtagccatg  
ccctgttgaagatgatctactgtccactgcccgggtctcgtgactgtgaagccatgttacaactactgtcacaacatcatgagaggtctgtttggccaaccaagggtatctgatttgaatgggaacaatttcatagatgctatgctgatgtgtggcagagaggttagaggtctttcaacattg  
aatcgggtcatggtccatcagatgtgaagatttctgactattatgaacatgcaggataatagtgttcaagtgctcagaaggttttccaggatgtggaccccccaagccccctccagctggacgaatttctcgttccatctctgaaagtgccttcagtgctcgttcagaccatcacccccga  
ggaacgccccaccacagcagctggcactagtgttgaccgactgttactgtatgtcaaggagaactgaacagggccaagaatttctgctccctccgagcaacgtttgcaacgatgagaggtgtgtcaggaacggcaatgaggatgactgttgaatgggaaaggcaaaagc  
aggtacctgtttgacgtgacaggaatggattagccaaccaggggcaacaaccagaggtccaggttgacaccagcaaacacacatactatcttctgcaaatcatgtgctctcagtgatgaccagcaagatgaagaatgcatacaatgggaacgacgtggactctttgatcatgtga  
tgaagtagtggagaagggaagtggaggtgctgtgagtacagcagtgccctcagagtttgactacaatgccactgaccatgctgggaagagtgcgaatgagaaggccagacgtgtgtgtcctcgtgggacagggcctacctcctcactgtctctgcatctgttctgttatgca  
gagagatggagataa