

PSMB10-proteasome (prosome, macropain) subunit, beta type, 10 Gene

Reference: PTXBC017198

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	LMP10; MECL1; beta2i; proteasome subunit beta type-10; low molecular mass protein 10; macropain subunit MECL-1; multicatalytic endopeptidase complex subunit MECL-1; proteasome (prosome, macropain) subunit, beta type, 10; proteasome MECL-1; proteasome cata		

For research use only.

GENE DESCRIPTION	proteasome (prosome, macropain) subunit, beta type, 10
GENE ID	5699
ORIGIN SPECIES	Human
ACCESSIONS	BC017198
VECTOR	pDONR223
NCBI SYMBOL	PSMB10
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgctgaagccagccctggagccccgaggggctctctctcagagaactgccaaagaaatgcatattggaacgcgtctcccggggctcaaggctccctacgcacgcaagaccgggaccaccatcgcgggctgtgttccaagacggggctattctgggcgccgatacgcgagcc
actaacgattcggtcgtggcgacaagagctgcgagaagatccacttcacgccccaaaatctactgctgtggggctggagtagccgcggagccgagatgaccacacggatgggtggcgctccaagatggagctacacgcgctatctacggccgcgagccccgctggccacggtc
actcgcatctgcgcagacgctcttcaggctaccagggccacgtgggtgcatcgctgatcgtggcgccgtgtagacctgactggaccgcagctctacgggtgtgcatcccatggctctacagccgtctgcccctacagccctgggctctggtcaggacgcggccctggcggtgtaga
agaccgggtccagccgaacatgacgctgagggctgctcaggggctgctggtggaagccgtcaccgccgggattctgggtgacctgggctccgggggcaatgtggacgcatgtgtatcacaaagactggcgccaagctgtcgcgacactgagctacccacagagccgtgaaga
ggctctggccgctaccactttgtgctggaaccacagctgctgacccagacagtgaagccactaacctggagctagtggaggaaactgtgcaggctatggaggtgagtaa