

C1orf128-chromosome 1 open reading frame 128 Gene

Reference: PTXBC017208

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C1orf128; HT014; TXNL1CL; PITH domain-containing protein 1; PITH (C-terminal proteasome-interacting domain of thioredoxin-like) domain containing 1; TXNL1 C-terminal like; PITH domain containing 1		

For research use only.

GENE DESCRIPTION	chromosome 1 open reading frame 128
GENE ID	57095
ORIGIN SPECIES	Human
ACCESSIONS	BC017208
VECTOR	pDONR223
NCBI SYMBOL	C1orf128
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtcgcacggtcacagccacggcggggtggctgccgtgcgccgaacgggaggagccgcccagcagcggcgccctggcctacggcctgtacctgcgcacgacctggagcggctgcaatgccttaacgagagccgcgaggcagcggcgccggtcttcaagccgtg
ggaggagcggaccgaccgctccaagttgttgaagtgatgcagatgaagagctctgtttaatatctcatttacggcgaatgtcaagctcaaggcatcattataatggagaggatgatgactcacaccctctgagatgagactgtacaagaatattccacagatgtccttgatgatacag
aaaggagccagatcagacctttagctgaaccgggatcttacaggagaattagagatgctacaaaaattctcgtttttcaaatgtctatcatctcaattcatattcaaaaaacttcggagcagatagacaaaaggcttttatattggcctgagaggagagtgaggactgagcttcgccgacac
gaggtgacctctgcaattacgaagcatctgccaaccagcagaccatagggtccatcagggttaccaccacagacacctttattcttaa