



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

RNF166-ring finger protein 166 Gene

Reference: PTXBC017226

Product type:

DNA

Storage Condition:

-20°C

Size/Quantity:

2 ug

Delivery Condition:

+4°C, Blue Ice

Synonyms

RING finger protein 166

For research use only.

GENE DESCRIPTION	ring finger protein 166
GENE ID	115992
ORIGIN SPECIES	Human
ACCESSIONS	BC017226
VECTOR	pDONR223
NCBI SYMBOL	RNF166
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgctgatcttcgcagcctggtggtccctcgctcagcagcggcagcgcgcggcgccggcgccggcgccgcgacagcgctcgaggcgcagtaacctgcccatctccttgagggtctatcaccggcccggtgccatcgcgcgtcgccacacgttctcggggagtgctctcc
agccctcgctcaggtgtccatccccgtctgtcccactgtgccgctgcttctgaccccaagaaggtgtggaacaggccaccacgttgagaaagcagctctcatctcaaaagccctgtcgcaggctcgaacaaaaaggtgaccttgcgaagatgagatgcacatttgcctctcgt
aaggttcaggagcagatgtgcccaactgcccaagtctgtccctcgtgtgtgcccaatcacagcctatcccagcaacatcccccaaggtcaccttgcctgcctgtatctgtgtgctccgcaactcgtgaccacagaggactgtgtgaagcactgtgtgtgaagaccacgcgacgacccc
aacccgtgtgtgtgccccatctcgtcgcgaatgccctggggggagccccagctacaagagcgcgaacttctcgcagcactgtctaccgcacaagttctctacgacaccttgttgactacagtattgacgaggaggccgccttcaggctcgtctgtgccctgtctctctgagaactg
a

a

More info: Please contact us for any question: contact@proteogenix.fr

Page 1/1