



Reference: PTXBC017238

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	tetatricopeptide repeat protein 17; TPR repeat protein 17; tetatricopeptide repeat domain 17		

For research use only.

GENE DESCRIPTION	tetratricopeptide repeat domain 17
GENE ID	55761
ORIGIN SPECIES	Human
ACCESSIONS	BC017238
VECTOR	pDONR223
NCBI SYMBOL	TTC17
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]