



Reference: PTXBC017256

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HKR1, GLI-Kruppel zinc finger family member; oncogene HKR1; GLI-Kruppel family member HKR1; ZNF875; Krueppel-related zinc finger protein 1; zinc finger protein 875		

For research use only.

GENE DESCRIPTION	GLI-Kruppel family member HKR1
GENE ID	284459
ORIGIN SPECIES	Human
ACCESSIONS	BC017256
VECTOR	pDONR223
NCBI SYMBOL	HKR1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]