

TLX3-T-cell leukemia homeobox 3 Gene

Reference: PTXBC017291

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	RNX; T-cell leukemia homeobox protein 3; homeo box 11-like 2; homeobox protein Hox-11L2; T-cell leukemia homeobox 3		

For research use only.

GENE DESCRIPTION	T-cell leukemia homeobox 3
GENE ID	30012
ORIGIN SPECIES	Human
ACCESSIONS	BC017291
VECTOR	pDONR223
NCBI SYMBOL	TLX3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaggcgcccgccagcgcagaccccgaccgcagcccatcagctcggcatcaccagatccttaacagcccgaccagcagcgcacccgcccccgggggccccgacggcgccagctacctgggagggcccccgggggccgtccggcgccacataccg
tctctgcccgcctctttgcgccctcggcgcccttcgaggacgaggatcttacagtgtaacctgagcctagcggcgcccgaggcggtgatccgggtgcccgcgcacaggccgctgccccggggccgtgccaccgcctctgccaagcgcgctaccgcccattgccctccgtgccacg
gtctccagccttgccggtctcaatttccccctggatggagagcagccgcttctgtaaaagaccgcttcacagcggcgccgcactcagcccttcaccgtgacccggcgcatcgccacccctaccagaaccggagcgcgccaagcgttaagaagccgcacgtctttccggg
tgcagatctcgagctgaaaaagcgcttcctcgcagagaagctacgtgcccctctgcccagagggcgccgctcgccaagtccctcaaaatgacggacgcgcaggtcaagacgtggttccaaaacggagaccgaagtggcgccgagacggcgaggagcgggagggcgagcgg
cagcaggcgagccggtcatgctgagctgcaaacagcgccttccaaaagagcctcaacgactccatccagcctgacccgctctgtctgcacaactcgtcactcttctgtctgcagaatctgcagccctgggagggagatgtccaaggttcccgtgtcacctccctggtgtga