

ZNHIT1-zinc finger, HIT type 1 Gene

Reference: PTXBC017333

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CG1I; ZNFN4A1; zinc finger HIT domain-containing protein 1; H_DJ0747G18.14; cyclin-G1-binding protein 1; p18 Hamlet; p18Hamlet; zinc finger protein subfamily 4A member 1; zinc finger protein, subfamily 4A (HIT domain containing), member 1; zinc finger, HI		

For research use only.

GENE DESCRIPTION	zinc finger, HIT type 1
GENE ID	10467
ORIGIN SPECIES	Human
ACCESSIONS	BC017333
VECTOR	pDONR223
NCBI SYMBOL	ZNHIT1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggtggagaagaacttcggttcgctccaggaccccggcagcggcggtgctggaccgggctgccggcagcgtcgcataaccggcagctggaggccctggagaatgacaacttcaggatgacccccacggggactccctcagctcggcaagagactgcctcagttga
tgacgafcgacactggaagaaaaaagaaaaacccgaggtgatcattttaacttcgcttccgaaaaaacttcaggccctgttggaggagcagaacttgagtgtgcccaggggccctaactacctgacggcctgtgcgggacccccatcgcccccagcgccttctgtgctgt
ctgtggttccatccccctacacctgtgtcagctgcgggtgccggtactgcactgtgcgctgtctggggaccaccaggagaccaggtgtctgaagtggactgtgtga