

RPL26L1-ribosomal protein L26-like 1 Gene

Reference: PTXBC017360

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	RPL26P1; 60S ribosomal protein L26-like 1; ribosomal protein L26 homolog; ribosomal protein L26 pseudogene 1; ribosomal protein L26 like 1		

For research use only.

GENE DESCRIPTION	ribosomal protein L26-like 1
GENE ID	51121
ORIGIN SPECIES	Human
ACCESSIONS	BC017360
VECTOR	pDONR223
NCBI SYMBOL	RPL26L1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaagttcaatcccttctgttacctcggaccgcagtaaaaaacgcaaacgtcaattcaatgccccctcacacgtgcgcaggagatcatgtcatcccgctctccaaggagctgcgcagaagtacaatgtccgctccatgccatccgcaaggacgacgaggtccaggtagtcgagga
cactacaaagggtcagcaaatggcaaggtagtccaggtgtacagaaagaatatgtcatctacatcgagcgggtgcagcgtgagaaggccaacggcacaactgtccacgtgggcattcaccacccaagcaggtgggtatccaggctaaaactggacaaggatcggaaaaaattcttga
acgcaaagccaagtcgcagaagttgaaaaagagaaggcaaatataagaagaacttattgagaaaatgcaggaataa