

LRP12-low density lipoprotein-related protein 12 Gene

Reference: PTXBC017381

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ST7; low-density lipoprotein receptor-related protein 12; suppressor of tumorigenicity 7 protein; LDL receptor related protein 12		

For research use only.

GENE DESCRIPTION	low density lipoprotein-related protein 12
GENE ID	29967
ORIGIN SPECIES	Human
ACCESSIONS	BC017381
VECTOR	pDONR223
NCBI SYMBOL	LRP12
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaagatattctactatgcattaaatctgaactttataaaacatgtacaaaaattgtacaagataagttccacctgtaattgtcttccctaaacagggttcgcttcattggaccctagggttgcactaaaattatatacaagggtcagatgagcttagtcacaaagcactatcactttaaac
tattatttgcaccacagcaactatataattccatagcttttgctggggcgggggacattttattacaactgaaattgcttgcgtgtttcatatttattgtgtatttaaaaaatacattgtgttaagagtgatttttcaatataattttattctggggggatcatgctacactctcaaaagaaattaa