

NT5C-5', 3'-nucleotidase, cytosolic Gene

Reference: PTXBC017454

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	DNT; DNT1; HEL74; P5N2; PN-I; PN-II; UMPH2; cdN; dNT-1; 5'(3')-deoxyribonucleotidase, cytosolic type; 5' nucleotidase, deoxy (pyrimidine), cytosolic type C; cytosolic 5',3'-pyrimidine nucleotidase; deoxy-5'-nucleotidase 1; epididymis luminal protein 74; u		

For research use only.

GENE DESCRIPTION	5', 3'-nucleotidase, cytosolic
GENE ID	30833
ORIGIN SPECIES	Human
ACCESSIONS	BC017454
VECTOR	pDONR223
NCBI SYMBOL	NT5C
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgcggagcgtgcgcgtgctgggacatggacggcgtctggccgacttcgagggcggcctcctgcggggcttcgcgcgcgttcctgaggagccgcacgtgccgctggagcagcgcgcggttcctggcccgagcagtagccgcctgcggcccgacctggcg
gataaagtggccagtggtgacgaagcccccggccttttctggacctggagccatccccgggagccttggacgtgtgcgggagatgaacgacctaccggacacgcaggtcttcatctgcaccagccccctgctgaagtaccaccactgtgtgggtgagaagtaccgtgggtggagca
gcacctggggcccgatcgtagaacgaattatcctgacaagggaagacgggtgtcttgggggacctgctcattgatgacaaggacacagttcgagggcagggagacccaagctgggagcacatcttcttcacctgtgccacaatcggcacctgggtcctgccccgacaagg
agacggctgctcctggagtacaactggaggagatcttagatagcaagcgcggagctgcgcagcgggaatga