

C7orf45-chromosome 7 open reading frame 45 Gene

Reference: PTXBC017587

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C7orf45; serine-rich single-pass membrane protein 1; serine rich single-pass membrane protein 1		

For research use only.

GENE DESCRIPTION	chromosome 7 open reading frame 45
GENE ID	136263
ORIGIN SPECIES	Human
ACCESSIONS	BC017587
VECTOR	pDONR223
NCBI SYMBOL	C7orf45
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggagacctttttcttattttgggaggtagatcctccccaatacctgttaattgtgccattccaaatcaggattatgaatgctggaagatgactctgtggaaccataggagcttctgctttggtattttgtatctatttgcctgatgttcttcttagggcttctgtctggtatgtctgagga
taaaaaggatgaaggcagtgaggacaagtactcagtaaggaaagcaagcaagagacttctgtgaagcgcaaaagacagtgccctgggatccctcacaacaatgaagaaccaaagcagaacaaactaccctgtaccaactcagaagtggtttgtcaatgcttatctg
aacaagacgagccaggcgccagtctcagttcaatgaggtgaaccagaaccaatgacagtgaatactacgagatgagcagtgaaagtctaaactcagaagcctctctggaaggagagtgaaagtgaacaccaccatcaccagacagtattaagaggagaaaaatggctcaga
ggcaagggaatctgggaagtacaaatgagcgaaggcactgctccactgcaagccttgagaaccaacgaatggtggcgccaccattccgacagaagccttcagtaacaccgccaatgaaaagagacagtcagaggaaagttccatctgacattacaagaatttagtaatt
ttga