



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

CXCL16-chemokine (C-X-C motif) ligand 16 Gene

Reference: PTXBC017588

Product type:

DNA

Storage Condition:

-20°C

Size/Quantity:

2 ug

Delivery Condition:

+4°C, Blue Ice

Synonyms

transmembrane chemokine CXCL16; CXCLG16; SR-PSOX; SRPSOX; C-X-C motif chemokine 16; CXC chemokine ligand 16; chemokine (C-X-C motif) ligand 16; scavenger receptor for phosphatidylserine and oxidized low density lipoprotein; small-inducible cytokine B16; C

For research use only.

GENE DESCRIPTION	chemokine (C-X-C motif) ligand 16
GENE ID	58191
ORIGIN SPECIES	Human
ACCESSIONS	BC017588
VECTOR	pDONR223
NCBI SYMBOL	CXCL16
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr

Page 1/1