

DNAJB3-DnaJ (Hsp40) homolog, subfamily B, member 3 Gene

Reference: PTXBC017590

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HCG3; dnaJ homolog subfamily B member 3; DnaJ (Hsp40) homolog, subfamily B, member 3; DnaJ heat shock protein family (Hsp40) member B3		

For research use only.

GENE DESCRIPTION	DnaJ (Hsp40) homolog, subfamily B, member 3
GENE ID	414061
ORIGIN SPECIES	Human
ACCESSIONS	BC017590
VECTOR	pDONR223
NCBI SYMBOL	DNAJB3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggtggactactacgagtgctggacgtgccccggcaggcctatccgagccatcaagaaggcgtaccgcaagctggcgctcaagtggcaccgcacaaaacccctgagaacaaggaggaagcggagaggagattcaagcagtgccgaggcctacgaggtgtgtcggac
gccaagaacgcgatatctatgaccgtatggcgagggcgcgaggcggtgcacaggcggcagggcccttcgaggacccttcgagttacgtcttcagcttcgcgaccagccgacgtcttcaggaggtcttcggcgccaggaccattctcttgacctctgggaaaccc
gctggagaatatttgggggggtcagaggaactgctggggaagcagaagcagagcgtctgcacccctttctctgccttcagtga