



## Reference: PTXBC017707

|                |                                                                                                                                                             |                     |                |
|----------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                                                                         | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                                                                        | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | FXY2; MRX101; RNF60; TRIM1; RING finger protein 60; midin 2; midline defect 2; tripartite motif protein 1; tripartite motif-containing protein 1; midline 2 |                     |                |

**For research use only.**

|                         |                                                                                                  |
|-------------------------|--------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | midline 2                                                                                        |
| <b>GENE ID</b>          | 11043                                                                                            |
| <b>ORIGIN SPECIES</b>   | Human                                                                                            |
| <b>ACCESSIONS</b>       | BC017707                                                                                         |
| <b>VECTOR</b>           | pDONR223                                                                                         |
| <b>NCBI SYMBOL</b>      | MID2                                                                                             |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                   |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

[illegible]

More info: Please contact us for any question: [contact@proteogenix.fr](mailto:contact@proteogenix.fr)