



Reference: PTXBC017715

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CSCF; FMD2; MEKK7; TGF1a; mitogen-activated protein kinase kinase kinase 7; TGF-beta activated kinase 1; transforming growth factor-beta-activated kinase 1		

For research use only.

GENE DESCRIPTION	mitogen-activated protein kinase kinase kinase 7
GENE ID	6885
ORIGIN SPECIES	Human
ACCESSIONS	BC017715
VECTOR	pDONR223
NCBI SYMBOL	MAP3K7
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr