



## Reference: PTXBC017727

|                               |          |                     |                |
|-------------------------------|----------|---------------------|----------------|
| Product type:                 | DNA      | Storage Condition:  | -20°C          |
| Size/Quantity:                | 2 ug     | Delivery Condition: | +4°C, Blue Ice |
| Synonyms                      | rhoteKin |                     |                |
| <b>For research use only.</b> |          |                     |                |

|                         |                                                                                                   |
|-------------------------|---------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | rhotekin                                                                                          |
| <b>GENE ID</b>          | 6242                                                                                              |
| <b>ORIGIN SPECIES</b>   | Human                                                                                             |
| <b>ACCESSIONS</b>       | BC017727                                                                                          |
| <b>VECTOR</b>           | pDONR223                                                                                          |
| <b>NCBI SYMBOL</b>      | RTKN                                                                                              |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                    |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

**ORF (Only for reference, product provided without STOP codon )**

[illegible]

More info: Please contact us for any question: [contact@proteogenix.fr](mailto:contact@proteogenix.fr)