



Reference: PTXBC017758

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CYPIVB1; P-450HP; cytochrome P450 4B1; cytochrome P450, family 4, subfamily B, polypeptide 1; cytochrome P450, subfamily IVB, polypeptide 1; cytochrome P450-HP; microsomal monooxygenase; cytochrome P450 family 4 subfamily B member 1		

For research use only.

GENE DESCRIPTION	cytochrome P450, family 4, subfamily B, polypeptide 1
GENE ID	1580
ORIGIN SPECIES	Human
ACCESSIONS	BC017758
VECTOR	pDONR223
NCBI SYMBOL	CYP4B1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr