

C3orf14-chromosome 3 open reading frame 14 Gene

Reference: PTXBC017772

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	uncharacterized protein C3orf14; HT021; chromosome 3 open reading frame 14		
For research use only.			

GENE DESCRIPTION	chromosome 3 open reading frame 14
GENE ID	57415
ORIGIN SPECIES	Human
ACCESSIONS	BC017772
VECTOR	pDONR223
NCBI SYMBOL	C3orf14
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgacttcctgtttgctcaagaaattcgctttctaaaagacatgaagaaatagatcacaaagattaatgttacttcaacaatggagaataaattgggtgatcaacacagaaaaggcatctcaactccaactgttgagactgctttaaaaggaaaccttagtcttttaaggatatagaagcagcagaaaagtcactacagaccaggattcaccacattccacggcctgaggtgtttctcttgagactcgttactggcatcagtagaagaatatattccaaatgggaacagtgttcttttaggaagagaccataatcctttgctgttgaaaatcaaatgaagcagaaaataccattcaaatgaggcacagcgataa