

MPZL2-myelin protein zero-like 2 Gene

Reference: PTXBC017774

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	EVA1; myelin protein zero-like protein 2; epithelial V-like antigen 1; myelin protein zero like 2		

For research use only.

GENE DESCRIPTION	myelin protein zero-like 2
GENE ID	10205
ORIGIN SPECIES	Human
ACCESSIONS	BC017774
VECTOR	pDONR223
NCBI SYMBOL	MPZL2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtatggcaagagctctactctgtcgggtgtctctctctggcctacagctcacagctcttggcctatagcagctgtggaaattatctcccggtgtcggagctgttaatgggacagatgctcggtaaaatgcactttccagcttggccctgtgggtgatgcttaacagtgacctg
gaattttctctctagacgggggacctgagcagttgtatctactaccacatagatccctccaacccatgagtgggcgggttaaggaccgggtgtcttgggatgggaatcctgagcgggtacgatgctccatcctctctggaaactgcagttcgacgacaatgggacatacactgccag
gtgaagaaccaccctgatgttgatgggggataggggagatccggctcagcgtcgtgcacactgtacgtctctctgagatccacttctgctctggccattggctctgcctgtgactgatgataatagtaattgtagtgtctctccagcattaccgaaaaagcgaatggggcga
gagctcataaagtgggtggagataaaatcaaaagaaggaaagctcaaccaagagaaaaagggtctctgtttatttagaagacacagactaa