

CYTH4-cytohesin 4 Gene

Reference: PTXBC017780

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CYT4; DJ63G5.1; PSCD4; cytohesin-4; PH, SEC7 and coiled-coil domain-containing protein 4; pleckstrin homology, Sec7 and coiled/coil domains 4		

For research use only.

GENE DESCRIPTION	cytohesin 4
GENE ID	27128
ORIGIN SPECIES	Human
ACCESSIONS	BC017780
VECTOR	pDONR223
NCBI SYMBOL	CYTH4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggacctgtgccaccagagccccggagctgagcagcggggagacggaaggttacagagatcaagtggcaccgaaagcagctcctggagacatccagaagctgaaggatgagattgcagatgtttgcccacatcgactcttcgagagtcggaggagagccggatgg
cccagaaggagaaggagctgtgtattggcgcaagaagtcaacatggaccccgccaagggtatccagttatctatgagcacaagctgctgacccctgacgtccaggacattgcacgggttcctgtataaaggcgaggcgctcaacaagacagccattggtacctacctggggagagg
taa