



Reference: PTXBC017782

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CCN5; CT58; CTGF-L; WNT1-inducible-signaling pathway protein 2; CCN family member 5; connective tissue growth factor-like protein; connective tissue growth factor-related protein 58; WNT1 inducible signaling pathway protein 2		

For research use only.

GENE DESCRIPTION	WNT1 inducible signaling pathway protein 2
GENE ID	8839
ORIGIN SPECIES	Human
ACCESSIONS	BC017782
VECTOR	pDONR223
NCBI SYMBOL	WISP2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]